

tün Ultrasonics

Original tün® Ultrasonic Tips (E1-E6) are made from all stainless steel, providing the strength and durability that our clinicians have come to rely on with every procedure. Their unique design allows for greater radial tip movement and control than conventional ultrasonic tips.

Our growing variety of tün® ultrasonics featuring tün Gold B, C, K lines, as well as our FRR lines offer the exact dimensions of competitors' tips at a fraction of the cost.

7mm of diamond coating, hardened by a special Duracoat™abrasive empowers tün® with exceptional cutting performance.

All tün® tips come in both 3mm x 0.6mm (Satelec thread style) as well as 3mm x 0.5mm (EMS thread style), and fit a plethora of endodontic ultrasonic units from all the leading manufacturers, including:

- Satelec Units
- NSK
- ADEC
- Woodpecker
- Obtura Spartan
- Dentsply/Tulsa
- ASI
- J. Morita
- Sybron Endo
- Vista

tün Ultrasonics

by **engineered
endodontics®**

Call:
(262) 501-0075

or visit:
www.EngineeredEndo.com

tün Ultrasonics

ActaSonic

File Retrieval

Endodontic

Restorative

Gold B/C/K



ActaSonic

TUN-AS ActaSonic Irrigation Tip

ActaSonic Irrigation Tip is a solid single piece titanium tip for easy navigation of the canal, and provides maximum tip movement for optimized agitation of canal debris and solution movement without removing dentin or enlarging the canal.

File Retrieval (FRR/FRR-D)

FRR 8 27mm Titanium

FRR 7 24mm Titanium

FRR 6 20mm Titanium

Titanium FRR tips have the flexibility & strength of titanium alloy to remove obstructions in the canal, separated instruments and aid in retreatments.

FRR 8-D 27mm Diamond-coated Titanium

FRR 7-D 24mm Diamond-coated Titanium

FRR 6-D 20mm Diamond-coated Titanium

Tün FRR-D has all the abilities of titanium alloy tips, but with micro diamond coating for a smoother, more efficient cutting performance.

Endodontic

E1 - Taper Tip w/Diamond

E3 - Ball Tip

E4 - Football Tip

E5 - Cone Tip

E2 - Taper Tip

E6 - Post Buster

Tün® Endodontic tips are made from all stainless steel, armored with Duracoat diamond plating for durable, efficient cutting performance, and are available in Taper (with & without diamond), Ball, Cone, Football and Post Buster styles.

Restorative

R3 (1.1mm)

R2 (1.3mm x 1.7mm taper)

R1 (1mm x 1.4mm taper)

Tün® Restorative Tips are crafted for an accurate pre-prosthetics finish, and precise preparation for crowns.

Gold B/C/K

GB-1

GB-1A

GB-2

GB-2A

GB-3

GB-3A

GC-5D (24 mm)

GC-4D (19 mm)

GC-3D (17 mm)

GC-2D (17 mm)

GC-1 (17 mm)

C-Line

Designed for retreatment with exact specs to CPR tips.

B-Line

Dimensionally identical to BUC tips.

K-Line

Dimensionally identical to KiS tips.

GK-2D GK-1D GK-6D GK-5D GK-3D GK-4D